



**Maintenance Solutions**

## **SUMMARY OF ANTIMICROBIAL ACTIVITY**

### **A2Z DISINFECTING GLASS & MULTISURFACE CLEANER**

#### **Description**

A2Z is a one-step Ready-To-Use hospital disinfectant cleaner that is effective against a broad spectrum of bacteria, viruses or fungi and will inhibit the growth of mold and mildew and their odors. It uses 5th generation quaternary ammonium chlorides as its disinfectant base that has been proven highly effective in the presence of organic soils.

---

#### **Regulatory Summary**

<b>EPA Registration No.</b>	10324-85-8325
<b>USDA Authorization</b>	None
<b>California Status</b>	None
<b>Canadian PCP#</b>	None
<b>Canadian Din #</b>	None

#### **Physical Properties**

<b>pH as is</b>	7.0 - 9.0	<b>Flash Point (PMCC)</b>	None
<b>Specific Gravity @ 25°C</b>	1.000	<b>% Quat (mol. wt.360.5)</b>	0.086 - 0.0946
<b>Pounds per gallon @25°C</b>	8.36	<b>% Volatile</b>	98+

## Efficacy

### Hospital Disinfection

This product is bactericidal according to the AOAC Use Dilution Test method on hard inanimate surfaces modified in the presence of 5% organic serum (850 ppm active). Treated surfaces must remain wet for 10 minutes

(Testing is performed per the AOAC UDT/GST method (DIS/TSS-1). Sixty carriers are required on 3 separate lots, one of which must be > 60 days old against *Pseudomonas aeruginosa*, *Salmonella enterica* and *Staphylococcus aureus*. Killing of 59 out of 60 carriers is required (total carriers = 540).)

Organism	Carrier Population	Sample	# Carriers	# Positive
<i>Pseudomonas aeruginosa</i> ATCC #15442	2.3 X 10 <sup>6</sup> CFU/Carrier	A (60 Days Old)	60	0/60
	1.7 X 10 <sup>6</sup> CFU/Carrier	B	60	0/60
	1.3 X 10 <sup>6</sup> CFU/Carrier	C	60	1/60
<i>Salmonella enterica</i> ATCC #10708	1.1 X 10 <sup>5</sup> CFU/Carrier	A (60 Days Old)	60	0/60
	1.5 X 10 <sup>6</sup> CFU/Carrier	B	60	0/60
	2.1 X 10 <sup>6</sup> CFU/Carrier	C	60	0/60
<i>Staphylococcus aureus</i> ATCC #6538	1.5 X 10 <sup>6</sup> CFU/Carrier	A (60 Days Old)	60	0/60
	1.4 X 10 <sup>6</sup> CFU/Carrier	B	60	0/60
	4.7 X 10 <sup>5</sup> CFU/Carrier	C	60	0/60

### Supplemental Organisms

(Testing is performed per the AOAC UDT/GST method. Ten carriers are required on 2 separate lots against each supplemental organism. Killing of 10 out of 10 carriers is required (total carriers = 20).)

Organism	Carrier Population	Sample	# Carriers	# Positive
<i>Burkholderia cepacia</i> ATCC 25416	3.5 X 10 <sup>6</sup> CFU/Carrier	A	10	0/10
		B	10	0/10
<i>Campylobacter jejuni</i> ATCC 29428	2.9 X 10 <sup>5</sup> CFU/Carrier	A	10	0/10
		B	10	0/10
<i>Corynebacterium ammoniagenes</i> ATCC 6871	1.8 X 10 <sup>5</sup> CFU/Carrier	A	10	0/10
		B	10	0/10
<i>Enterobacter aerogenes</i> ATCC 13048	4.1 X 10 <sup>6</sup> CFU/Carrier	A	10	0/10
		B	10	0/10
<i>Enterobacter cloacae</i> Clinical Isolate	3.9 X 10 <sup>5</sup> CFU/Carrier	A	10	0/10
		B	10	0/10
<i>Enterobacteriaceae</i> with extended beta-lactamase resistance ATCC BAA-72	6.6 X 10 <sup>4</sup> CFU/Carrier	A	10	0/10
	1.25 X 10 <sup>6</sup> CFU/Carrier	B	10	0/10
<i>Enterococcus faecalis</i> ATCC 19433	9.4 X 10 <sup>4</sup> CFU/Carrier	A	10	0/10
		B	10	0/10
<i>Enterococcus faecium</i> Vancomycin Resistant (VRE)	4.5 X 10 <sup>5</sup> CFU/Carrier	A	10	0/10
		B	10	0/10
<i>Escherichia coli</i> ATCC 11229	3.2 X 10 <sup>5</sup> CFU/Carrier	A	10	0/10
		B	10	0/10
<i>Escherichia coli</i> Antibiotic Resistant Clinical Isolate	3.9 X 10 <sup>5</sup> CFU/Carrier	A	20	0/20
		B	20	0/20

<i>Escherichia coli</i> 0157:H7 ATCC 35150	1.1 X 10 <sup>4</sup> CFU/Carrier	A	20	0/20
		B	20	0/20
<i>Klebsiella pneumoniae</i> ATCC 4352	9.9 X 10 <sup>4</sup> CFU/Carrier	A	10	0/10
		B	10	0/10
<i>Klebsiella pneumoniae</i> Antibiotic Resistant Clinical Isolate	2.7 X 10 <sup>5</sup> CFU/Carrier	A	10	0/10
		B	10	0/10
<i>Legionella pneumophila</i> ATCC 33153	8.2 X 10 <sup>7</sup> CFU/Carrier	A	10	0/10
		B	10	0/10
<i>Listeria monocytogenes</i> ATCC 984	1.85 X 10 <sup>5</sup> CFU/Carrier	A	10	0/10
		B	10	0/10
<i>Proteus mirabilis</i> Clinical Isolate	1.9 X 10 <sup>6</sup> CFU/Carrier	A	20	0/20
		B	20	0/20
<i>Proteus vulgaris</i> ATCC 33420	4.55 X 10 <sup>4</sup> CFU/Carrier	A	20	0/20
		B	20	0/20
<i>Pseudomonas aeruginosa</i> Clinical Isolate	1.2 X 10 <sup>6</sup> CFU/Carrier	A	10	0/10
		B	10	0/10
<i>Salmonella typhi</i> ATCC 6539	5.1 X 10 <sup>5</sup> CFU/Carrier	A	10	0/10
		B	10	0/10
<i>Serratia marcescens</i> ATCC 43861	1.5 X 10 <sup>5</sup> CFU/Carrier	A	10	0/10
		B	10	0/10
<i>Shigella dysenteriae</i> ATCC 9361	5.45 X 10 <sup>4</sup> CFU/Carrier	A	10	0/10
		B	10	0/10
<i>Shigella flexneri</i> ATCC 12022	4.85 X 10 <sup>4</sup> CFU/Carrier	A	20	0/20
		B	20	0/20
<i>Shigella sonnei</i> ATCC 9290	2.75 X 10 <sup>4</sup> CFU/Carrier	A	20	0/20
		B	20	0/20
<i>Staphylococcus aureus</i> (Methicillin Resistant) (MRSA) ATCC 33591	1.45 X 10 <sup>5</sup> CFU/Carrier	A	10	0/10
		B	10	0/10
Community Associates Methicillin Resistant <i>Staphylococcus aureus</i> (CA-MRSA) (NRS) (Genotype USA400)	2.77 X 10 <sup>5</sup> CFU/Carrier	A	10	0/10
		B	10	0/10
<i>Staphylococcus epidermidis</i> Antibiotic Resistant Clinical Isolate	4.2 X 10 <sup>5</sup> CFU/Carrier	A	10	0/10
		B	10	0/10
<i>Streptococcus pyogenes</i> ATCC 19615	3.35 X 10 <sup>6</sup> CFU/Carrier	A	10	0/10
		B	10	0/10
<i>Vibrio cholera</i> ATCC 11623	9.3 X 10 <sup>6</sup> CFU/Carrier	A	10	0/10
		B	10	0/10

Nothing contained herein grants or extends a license, express or implied, in connection with patents, issued or pending, of the manufacturer or others. The information contained herein is based on the manufacturer's own study and the works of others. The manufacturer makes no warranties, expressed or implied, as to the accuracy, completeness, or adequacy of the information contained herein. The manufacturer shall not be liable (regardless of fault) to the vendee's employees, or anyone for any direct, special or consequential damages arising out of or in connection with the accuracy, completeness, adequacy or furnishing of such information.

### Virucidal against

This product was evaluated in the presence of 5% serum (850 ppm quat active), with a 10 minute contact time and found to be effective against the following viruses on hard nonporous environmental surfaces.

Organism	Dried Virus Control;	Sample	Result	Log Reduction
Avian Influenza A (H5N1) Virus	4.5 Log <sub>10</sub>	A	≤0.5 Log <sub>10</sub>	≥4.0 Log <sub>10</sub>
		B	≤0.5 Log <sub>10</sub>	≥4.0 Log <sub>10</sub>
Avian influenza /Turkey/Wisconsin ATCC VR-798	6.0 Log <sub>10</sub>	A	≤1.5 Log <sub>10</sub>	≥4.5 Log <sub>10</sub>
		B	≤1.5 Log <sub>10</sub>	≥4.5 Log <sub>10</sub>
Canine Coronavirus ATCC VR-809	4.75 Log <sub>10</sub>	A	≤1.5 Log <sub>10</sub>	≥3.25 Log <sub>10</sub>
		B	≤0.5 Log <sub>10</sub>	≥4.25 Log <sub>10</sub>
Canine Distemper	5.0 Log <sub>10</sub>	A	≤1.5 Log <sub>10</sub>	≥3.5 Log <sub>10</sub>
		B	≤1.5 Log <sub>10</sub>	≥3.5 Log <sub>10</sub>
Hantavirus	5.0 Log <sub>10</sub>	A	≤1.5 Log <sub>10</sub>	≥3.5 Log <sub>10</sub>
		B	≤1.5 Log <sub>10</sub>	≥3.5 Log <sub>10</sub>
Hepatitis B Virus	5.5 Log <sub>10</sub>	A	≤1.5 Log <sub>10</sub>	≥4.0 Log <sub>10</sub>
	5.5 Log <sub>10</sub>	B	≤1.5 Log <sub>10</sub>	≥4.5 Log <sub>10</sub>
	4.5 Log <sub>10</sub>	Confirmatory A	≤1.5 Log <sub>10</sub>	≥3.0 Log <sub>10</sub>
Hepatitis C Virus ATCC CCL-22	6.84 Log <sub>10</sub>	A	≤1.51 Log <sub>10</sub>	≥5.33 Log <sub>10</sub>
	6.84 Log <sub>10</sub>	B	≤1.51 Log <sub>10</sub>	≥5.33 Log <sub>10</sub>
	7.14 Log <sub>10</sub>	Confirmatory B	≤1.7 Log <sub>10</sub>	≥5.44 Log <sub>10</sub>
Herpes Simplex Type 1 ATCC VR-260	5.0 Log <sub>10</sub>	A	≤1.5 Log <sub>10</sub>	≥3.5 Log <sub>10</sub>
		B	≤1.5 Log <sub>10</sub>	≥3.5 Log <sub>10</sub>
Herpes Simplex Type 2 ATCC VR-734	5.0 Log <sub>10</sub>	A	≤1.5 Log <sub>10</sub>	≥3.5 Log <sub>10</sub>
		B	≤1.5 Log <sub>10</sub>	≥3.5 Log <sub>10</sub>
Human Coronavirus ATCC VR-740	4.75 Log <sub>10</sub>	A	≤1.5 Log <sub>10</sub>	≥3.25 Log <sub>10</sub>
		B	≤1.5 Log <sub>10</sub>	≥3.25 Log <sub>10</sub>
Human Immunodeficiency Virus type 1 (HIV 1) HTLV-III <sub>RF</sub>	6.5 Log <sub>10</sub>	A	≤2.5 Log <sub>10</sub>	≥4.0 Log <sub>10</sub>
		B	≤2.5 Log <sub>10</sub>	≥4.0 Log <sub>10</sub>
<b>Influenza A (H1N1) virus ATCC VR-1469</b>	4.5 Log <sub>10</sub>	A	≤0.5 Log <sub>10</sub>	≥4.0 Log <sub>10</sub>
		B	≤0.5 Log <sub>10</sub>	≥4.0 Log <sub>10</sub>
Influenza A/Brazil Virus	4.8 Log <sub>10</sub>	A	≤0.5 Log <sub>10</sub>	≥4.3 Log <sub>10</sub>
		B	≤0.5 Log <sub>10</sub>	≥4.3 Log <sub>10</sub>

Infectious Bovine Rhinotracheitis virus (IBR) ATCC VR-188	5.0 Log <sub>10</sub>	A	≤1.5 Log <sub>10</sub>	≥3.5 Log <sub>10</sub>
		B	≤1.5 Log <sub>10</sub>	≥3.5 Log <sub>10</sub>
Newcastle disease virus ATCC VR-109	6.3 Log <sub>10</sub>	A	≤1.5 Log <sub>10</sub>	≥4.8 Log <sub>10</sub>
	5.8 Log <sub>10</sub>	B	≤1.5 Log <sub>10</sub>	≥4.3 Log <sub>10</sub>
Porcine Respiratory & Reproductive (PRRSV) Strain NVSL	5.5 Log <sub>10</sub>	A	≤1.5 Log <sub>10</sub>	≥4.0 Log <sub>10</sub>
		B	≤1.5 Log <sub>10</sub>	≥4.0 Log <sub>10</sub>
Porcine Rotavirus ATCC VR-893	4.5 Log <sub>10</sub>	A	≤1.5 Log <sub>10</sub>	≥3.0 Log <sub>10</sub>
		B	≤1.5 Log <sub>10</sub>	≥3.0 Log <sub>10</sub>
Pseudorabies virus ATCC VR-135	4.5 Log <sub>10</sub>	A	≤1.5 Log <sub>10</sub>	≥3.0 Log <sub>10</sub>
		B	≤1.5 Log <sub>10</sub>	≥3.0 Log <sub>10</sub>
Respiratory syncytial virus ATCC VR-26, Strain Long	4.5 Log <sub>10</sub>	A	≤1.5 Log <sub>10</sub>	≥3.0 Log <sub>10</sub>
		B	≤1.5 Log <sub>10</sub>	≥3.0 Log <sub>10</sub>
Transmissible Gastroenteritis (TGE) ATCC VR-742	5.7 Log <sub>10</sub>	A	≤2.5 Log <sub>10</sub>	≥3.2 Log <sub>10</sub>
		B	≤2.5 Log <sub>10</sub>	≥3.2 Log <sub>10</sub>
Vaccinia virus	5.5 Log <sub>10</sub>	A	≤1.5 Log <sub>10</sub>	≥4.0 Log <sub>10</sub>
		B	≤1.5 Log <sub>10</sub>	≥4.0 Log <sub>10</sub>

### Non-Food Contact Surface Sanitizer

This product is an effective on hard porous non-food contact surfaces. Treated surfaces must remain wet for 5 minutes. Then wipe with sponge, mop or cloth or allow to air dry. Food contact surfaces must be rinsed.

Testing is performed per EPA Guidance (DIS/TSS-10). Three lots are required, one of which must be ≥ 60 days old. Testing is performed against *Staphylococcus aureus* and *Klebsiella pneumoniae* containing 5% organic load. *Enterobacter aerogenes* may be substituted for *Klebsiella pneumoniae*. The results must show a reduction of at least 99.9% (3 Log<sub>10</sub>) in the number of each test microorganism over the parallel control count within 5 minutes.

Organism	Carrier Population	Sample	60 Second Kill cfu/Carrier	3 Minute Kill cfu/Carrier
<i>Klebsiella pneumoniae</i> ATCC 4352	5.46 Log <sub>10</sub>	A (60 Days Old)	3.69 Log <sub>10</sub>	4.55 Log <sub>10</sub>
		B	2.64 Log <sub>10</sub>	4.55 Log <sub>10</sub>
		C	4.46 Log <sub>10</sub>	4.55 Log <sub>10</sub>
<i>Staphylococcus aureus</i> ATCC #6538	5.18 Log <sub>10</sub>	A (60 Days Old)	3.84 Log <sub>10</sub>	5.26 Log <sub>10</sub>
		B	5.2 Log <sub>10</sub>	5.26 Log <sub>10</sub>
		C	5.2 Log <sub>10</sub>	5.26 Log <sub>10</sub>

### Mold and Mildew Control

Use this product to control the growth of mold and mildew and their odors on hard, non-porous surfaces. Thoroughly wet all treated surfaces completely. Let air-dry. Repeat application weekly or when growth or odor reappears.

Organism	Tile Number	Untreated After 7 Days	Sample A After 7 Days	Sample B After 7 Days
<i>Aspergillus niger</i> ATCC #6275	1	Growth 100%	No Growth 0%	No Growth 0%
	2	Growth 100%	No Growth 0%	No Growth 0%
	3	Growth 80%	No Growth 0%	No Growth 0%

	4	Growth 80%	No Growth 0%	No Growth 0%
	5	Growth 100%	No Growth 0%	No Growth 0%
	6	Growth 80%	No Growth 0%	No Growth 0%
	7	Growth 80%	No Growth 0%	No Growth 0%
	8	Growth 80%	No Growth 0%	No Growth 0%
	9	Growth 80%	No Growth 0%	No Growth 0%
	10	Growth 80%	No Growth 0%	No Growth 0%

### Fungicidal against

This product was evaluated at 2 ounces per gallon in the presence of 5% serum and 400 ppm hard water with a 10 minute contact time and found to be effective against the following fungi on hard nonporous environmental surfaces.

(Testing is performed per the AOAC fungicidal method (DIS/TSS-6). Two separate lots are tested against *Trichophyton mentagrophytes* in a suspension test. Killing of all fungal spores in 10 minutes is required.)

Organism	Carrier Population	Sample	# Carriers	# Positive
<i>Candida albicans</i> ATCC #10231	4.2 X 10 <sup>5</sup> CFU/Carrier	A	10	0/10
		B	10	0/10
<i>Dactylium dendroides</i> ATCC 6676	2.18 X 10 <sup>4</sup> CFU/Carrier	A	10	0/10
		B	10	0/10
<i>Trichophyton mentagrophytes</i> ATCC #9533	6.6 X 10 <sup>6</sup> CFU/Carrier	A	10	0/10
		B	10	0/10

The following data is for informational purposes only. This data was submitted to the EPA, but it was **not accepted** due to the human health issues that the agency has with the following pathogens. Even though this product was effective as shown the EPA **will not allow** these organisms to be added to the label.

Organism	Carrier Population	Sample	# Carriers	# Positive
<i>Penicillium variable</i> ATCC #32333	1.22 X 10 <sup>6</sup> CFU/Carrier	A	10	0/10
		B	10	0/10
<i>Stachbotrys chartarum</i> ATCC #66239	9.2 X 10 <sup>4</sup> CFU/Carrier	A	10	0/10
		B	10	0/10